

# EBP05054 - IMMUNIZING PEPTIDE

**SKU:** EBP05054



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_000256.4

**Immunogen** Peptide with sequence C-SESTKRKLASAV, from the C Terminus of the protein sequence according to NP\_000256.4.

**Product Comments** The peptide was used in the production of EB05054 - Goat Anti-NCF1 / p47phox Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-SESTKRKLASAV

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_plus\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png)

## DOCUMENTS

- [Data Sheet](#)