

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05067-B - IMMUNIZING PEPTIDE

SKU: EBP05067-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 002543.2; NP 001177810.1; NP 001177811.1 ID

Peptide with sequence C-DVRQWATSSLSWL., from the C Terminus of the protein sequence according to Immunogen

NP_002543.2; NP_001177810.1; NP_001177811.1.

The peptide was used in the production of EB05067-B - Goat Anti-ORC4L, Biotinylated Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-DVRQWATSSLSWL. Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade