

# EBP05067 - IMMUNIZING PEPTIDE

SKU: EBP05067



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_002543.2; NP\_001177810.1; NP\_001177811.1

**Immunogen** Peptide with sequence C-DVRQWATSSLSWL, from the C Terminus of the protein sequence according to NP\_002543.2; NP\_001177810.1; NP\_001177811.1.

**Product Comments** The peptide was used in the production of EBP05067 - Goat Anti-ORC4L Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-DVRQWATSSLSWL

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)