

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05073 - IMMUNIZING PEPTIDE

SKU: EBP05073

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_003141.2; NP_644805.1 ID

 $\begin{tabular}{l|l} \textbf{Immunogen} \\ \textbf{Peptide with sequence DMELTSECATSPM, from the C Terminus of the protein sequence according to NP_003141.2;} \\ \textbf{NP_644805.1.} \\ \end{tabular}$

The peptide was used in the production of EB05073 - Goat Anti-STAT3 (isoform 1 and 2) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

DMELTSECATSPM Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade