

EBP05074 - IMMUNIZING PEPTIDE

SKU: EBP05074



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060190.2; NP_001013863.1
Immunogen	Peptide with sequence C-LQKKLEKRRALEH, from the C Terminus of the protein sequence according to NP_060190.2; NP_001013863.1.
Product	The peptide was used in the production of EBP05074 - Goat Anti-STAP2 / BKS (aa391-403) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-LQKKLEKRRALEH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)