

EBP05079 - IMMUNIZING PEPTIDE

SKU: EBP05079



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_005301.2; NP_001025173.1

Immunogen Peptide with sequence LLRHCCTRVAL, from the C Terminus of the protein sequence according to NP_005301.2; NP_001025173.1.

Product Comments The peptide was used in the production of EBP05079 - Goat Anti-GRB7 (C Terminus) Antibody. It is currently untested in blocking applications.

Peptide Sequence LLRHCCTRVAL

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png