

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05080 - IMMUNIZING PEPTIDE

SKU: EBP05080

×

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

 $\textbf{Unit Size} \qquad 100 \; \mu \text{g}$

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_004481

 $\textbf{Immunogen} \ \ \textbf{Peptide with sequence KLKHYCARIAL, from the C Terminus of the protein sequence according to NP_004481.}$

Product The peptide was used in the production of EB05080 - Goat Anti-GRB14 Antibody. It is currently untested in

Comments blocking applications.

Peptide Sequence

KLKHYCARIAL

Shipping

Instructions Ambient

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$