

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05082 - IMMUNIZING PEPTIDE

SKU: EBP05082

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_004223 ID

Immunogen Peptide with sequence C-PNKMKDYLQEKHY, from the C Terminus of the protein sequence according to NP_004223.

Product The peptide was used in the production of EB05082 - Goat Anti-SOCS6 Antibody. It is currently untested in

Comments blocking applications.

Peptide

C-PNKMKDYLQEKHY Sequence

Shipping

Ambient Instructions

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png