

EBP05092 - IMMUNIZING PEPTIDE

SKU: EBP05092



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID AAF73914.1

Immunogen Peptide with sequence MPPAPSCPSERVTV, from the C Terminus of the protein sequence according to AAF73914.1.

Product Comments The peptide was used in the production of EBP05092 - Goat Anti-SH2B (683 aa isoform) Antibody. It is currently untested in blocking applications.

Peptide Sequence MPPAPSCPSERVTV

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png