



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05093 - IMMUNIZING PEPTIDE

SKU: EBP05093

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Immunogen Peptide with sequence C-GHPRAINNQYSFV, from the C Terminus of the protein sequence according to .

Product The peptide was used in the production of EB05093 - Goat Anti-SH2B (756 aa isoform) Antibody. It is currently

Comments untested in blocking applications.

Peptide

C-GHPRAINNQYSFV Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$ Grade