

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05094 - IMMUNIZING PEPTIDE

SKU: EBP05094

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_056318.2; NP_001295223.1 ID

Immunogen Peptide with sequence PFGASDCVTDHLP, from the C Terminus of the protein sequence according to NP_056318.2; NP_001295223.1.

The peptide was used in the production of EB05094 - Goat Anti-SH2-B / SH2B1 (isoform2) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

PFGASDCVTDHLP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade