

EBP05094 - IMMUNIZING PEPTIDE

SKU: EBP05094



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_056318.2; NP_001295223.1
Immunogen	Peptide with sequence PFGASDCVTDHLP, from the C Terminus of the protein sequence according to NP_056318.2; NP_001295223.1.
Product Comments	The peptide was used in the production of EBP05094 - Goat Anti-SH2-B / SH2B1 (isoform2) Antibody. It is currently untested in blocking applications.
Peptide Sequence	PFGASDCVTDHLP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiling_medium.png

DOCUMENTS

- [Data Sheet](#)