

EBP05100 - IMMUNIZING PEPTIDE

SKU: EBP05100



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_003020; NP_058544
Immunogen	Peptide with sequence SELCLQQPVERKQ, from the C Terminus of the protein sequence according to NP_003020; NP_058544.
Product	The peptide was used in the production of EB05100 - Goat Anti-SHC1; SHC3 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	SELCLQQPVERKQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png