

EBP05103 - IMMUNIZING PEPTIDE

SKU: EBP05103



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002342.1

Immunogen Peptide with sequence GIREDPDVCLKAP, from the C Terminus of the protein sequence according to NP_002342.1.

Product Comments The peptide was used in the production of EB05103 - Goat Anti-SH2D1A / SAP Antibody. It is currently untested in blocking applications.

Peptide Sequence GIREDPDVCLKAP

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png