

EBP05104 - IMMUNIZING PEPTIDE

SKU: EBP05104



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003868.1; NP_001257396.1
Immunogen	Peptide with sequence C-RLKDYLEEYKFQV, from the C Terminus of the protein sequence according to NP_003868.1; NP_001257396.1.
Product	The peptide was used in the production of EB05104 - Goat Anti-SOCS2 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-RLKDYLEEYKFQV
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png