

EBP05110 - IMMUNIZING PEPTIDE

SKU: EBP05110



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NM_007315

Immunogen Peptide with sequence C-GSVEFDSMMNTV, from the C Terminus of the protein sequence according to NM_007315.

Product The peptide was used in the production of EB05110 - Goat Anti-STAT1 alpha isoform Antibody. It is currently

Comments untested in blocking applications.

Peptide Sequence C-GSVEFDSMMNTV

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png