

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05111 - IMMUNIZING PEPTIDE**

**SKU:** EBP05111

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_037446.1; \ NP\_001107566.1; \ NP\_001245369.1; \ NP\_001245370.1; \ NP\_001245371.1$ ID

Peptide with sequence C-KDSTRLKYAVKVS, from the C Terminus of the protein sequence according to NP\_037446.1;

Immunogen NP\_001107566.1; NP\_001245369.1; NP\_001245370.1; NP\_001245371.1.

The peptide was used in the production of EB05111 - Goat Anti-BLNK / SLP-65 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

C-KDSTRLKYAVKVS Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade