

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05114 - IMMUNIZING PEPTIDE

SKU: EBP05114

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_002825.3; NP_001317366.1 ID

 $\label{eq:local_problem} \textbf{Immunogen} \ \ \ \ \ \ \ \ Peptide with sequence C-YENVGLMQQQKSFR, from the C Terminus of the protein sequence according to NP_002825.3; NP_001317366.1.$

Product The peptide was used in the production of EB05114 - Goat Anti-SHP2 / PTPN11 Antibody. It is currently untested in

Comments blocking applications.

Peptide

C-YENVGLMQQQKSFR Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade