

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05128 - IMMUNIZING PEPTIDE

SKU: EBP05128

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_071445 ID

Immunogen Peptide with sequence EEVDKLGCRDTRLLL, from the C Terminus of the protein sequence according to NP_071445.

Product The peptide was used in the production of EB05128 - Goat Anti-CARD15 / NOD2 (C Term) Antibody. It is currently

Comments untested in blocking applications.

Peptide

Sequence

EEVDKLGCRDTRLLL

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade