

EBP05138 - IMMUNIZING PEPTIDE

SKU: EBP05138



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_076926 and NP_057734
Immunogen	Peptide with sequence C-HMAETASPTPLNEV, from the C Terminus of the protein sequence according to NP_076926 and NP_057734.
Product	The peptide was used in the production of EB05138 - Goat Anti-MS4A4A Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-HMAETASPTPLNEV
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png