

EBP05139 - IMMUNIZING PEPTIDE

SKU: EBP05139



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NM_021201
Immunogen	Peptide with sequence CIQQVKSSSRSWI, from the C Terminus of the protein sequence according to NM_021201.
Product Comments	The peptide was used in the production of EB05139 - Goat Anti-MS4A7 (CFFM4) Antibody. It is currently untested in blocking applications.
Peptide Sequence	CIQQVKSSSRSWI
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png