

EBP05140 - IMMUNIZING PEPTIDE

SKU: EBP05140



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_068769.2; NP_690605.1
Immunogen	Peptide with sequence CQDQESSPIENDSSP, from the C Terminus of the protein sequence according to NP_068769.2; NP_690605.1.
Product Comments	The peptide was used in the production of EB05140 - Goat Anti-CD20 / MS4A1 (C Terminus) Antibody. It is currently untested in blocking applications.
Peptide Sequence	CQDQESSPIENDSSP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)