

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05140 - IMMUNIZING PEPTIDE**

**SKU:** EBP05140

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 068769.2; NP 690605.1 ID

Peptide with sequence CQDQESSPIENDSSP, from the C Terminus of the protein sequence according to Immunogen

NP 068769.2; NP 690605.1.

The peptide was used in the production of EB05140 - Goat Anti-CD20 / MS4A1 (C Terminus) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

CQDQESSPIENDSSP Sequence

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png