

EBP05150 - IMMUNIZING PEPTIDE

SKU: EBP05150



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Shipped at ambient temperature, store at -20°C.

Instructions

Immunogen Peptide with sequence PCWLEILLNNPR, from the C Terminus of the protein sequence according to .

Product The peptide was used in the production of EB05150 - Goat Anti-SMAD6 / MADH6 Antibody. It is currently

Comments untested in blocking applications.

Peptide
Sequence PCWLEILLNNPR

Shipping
Instructions Ambient

Product
Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png