

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05158 - IMMUNIZING PEPTIDE

SKU: EBP05158

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession $NP_057093.1; \ NP_001193580.1; \ NP_001193581.1; \ NP_001193582.1$ ID

Peptide with sequence C-QKGKVPITYLELLN, from the C Terminus of the protein sequence according to NP_057093.1;

Immunogen NP_001193580.1; NP_001193581.1; NP_001193582.1.

The peptide was used in the production of EB05158 - Goat Anti-BIF-1 / SH3GLB1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-QKGKVPITYLELLN Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade