

EBP05158 - IMMUNIZING PEPTIDE

SKU: EBP05158



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_057093.1; NP_001193580.1; NP_001193581.1; NP_001193582.1
Immunogen	Peptide with sequence C-QKGKVPITYLELLN, from the C Terminus of the protein sequence according to NP_057093.1; NP_001193580.1; NP_001193581.1; NP_001193582.1.
Product Comments	The peptide was used in the production of EB05158 - Goat Anti-BIF-1 / SH3GLB1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QKGKVPITYLELLN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)