

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05163 - IMMUNIZING PEPTIDE

SKU: EBP05163

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

AAA78935.1 ID

Peptide with sequence C-PGDVFYMPKE, from the internal region (near the C Terminus) of the protein sequence **Immunogen**

according to AAA78935.1.

The peptide was used in the production of EB05163 - Goat Anti-MTA1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-PGDVFYMPKE Sequence

Shipping Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade