

EBP05174 - IMMUNIZING PEPTIDE

SKU: EBP05174



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NM_021069

Immunogen Peptide with sequence SYYQRPFSAY-C, from the N Terminus of the protein sequence according to NM_021069.

Product Comments The peptide was used in the production of EB05174 - Goat Anti-ARGBP2 Antibody. It is currently untested in blocking applications.

Peptide Sequence SYYQRPFSAY-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png