

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05174 - IMMUNIZING PEPTIDE

SKU: EBP05174

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NM 021069 ID

Peptide with sequence SYYQRPFSPSAY-C, from the N Terminus of the protein sequence according to Immunogen

NM 021069.

The peptide was used in the production of EB05174 - Goat Anti-ARGBP2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

SYYQRPFSPSAY-C Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade