

EBP05186 - IMMUNIZING PEPTIDE

SKU: EBP05186



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004217.1

Immunogen Peptide with sequence SLPNPHELVSDLLC, from the C Terminus of the protein sequence according to NP_004217.1.

Product Comments The peptide was used in the production of EB05186 - Goat Anti-DRAK2 / STK17B Antibody. It is currently untested in blocking applications.

Peptide Sequence SLPNPHELVSDLLC

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png