

EBP05200 - IMMUNIZING PEPTIDE

SKU: EBP05200



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_115937.1; NP_742024.1
Immunogen	Peptide with sequence RALEDLLVSI GLDC, from the C Terminus of the protein sequence according to NP_115937.1; NP_742024.1.
Product	The peptide was used in the production of EB05200 - Goat Anti-BPOZ / ABTB1 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	RALEDLLVSI GLDC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)