

# EBP05218 - IMMUNIZING PEPTIDE

**SKU:** EBP05218



---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.                                                                                                                                                                  |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                            |
| <b>Storage</b>               |                                                                                                                                                                                                                                   |
| <b>Instructions</b>          | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                                   |
| <b>Accession ID</b>          | NP_001036194.1; NP_001036212.1; NP_005239.1                                                                                                                                                                                       |
| <b>Immunogen</b>             | Peptide with sequence C-TSAEPQYQPGDQT, from the C Terminus of the protein sequence according to NP_001036194.1; NP_001036212.1; NP_005239.1.                                                                                      |
| <b>Product</b>               | The peptide was used in the production of EB05218 - Goat Anti-FGR Antibody. It is currently untested in                                                                                                                           |
| <b>Comments</b>              | blocking applications.                                                                                                                                                                                                            |
| <b>Peptide Sequence</b>      | C-TSAEPQYQPGDQT                                                                                                                                                                                                                   |
| <b>Shipping Instructions</b> | Ambient                                                                                                                                                                                                                           |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)