



EBP05268 - IMMUNIZING PEPTIDE

SKU: EBP05268



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	AAF08298
Immunogen	Peptide with sequence DVRGSLQTPQNRPHC, from the internal region of the protein sequence according to AAF08298.
Product Comments	The peptide was used in the production of EB05268 - Goat Anti-SMURF1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	DVRGSLQTPQNRPHC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png