

EBP05273 - IMMUNIZING PEPTIDE

SKU: EBP05273



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_061846
Immunogen	Peptide with sequence C-QEARADLARRGLRF, from the C Terminus of the protein sequence according to NP_061846.
Product	The peptide was used in the production of EB05273 - Goat Anti-DPM3 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-QEARADLARRGLRF
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png