

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05278 - IMMUNIZING PEPTIDE

SKU: EBP05278

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_036269.1; NP_001127917.1; NP_001317584.1 ID

Immunogen Peptide with sequence TCCSVLINKKVDS, from the C Terminus of the protein sequence according to NP_036269.1; NP_001127917.1; NP_001317584.1.

The peptide was used in the production of EB05278 - Goat Anti-DDAH1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

TCCSVLINKKVDS Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade