

EBP05280 - IMMUNIZING PEPTIDE

SKU: EBP05280



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Shipped at ambient temperature, store at -20°C.

Instructions

Immunogen Peptide with sequence C-VKQAINKMLTKISS, from the C Terminus of the protein sequence according to .

Product The peptide was used in the production of EB05280 - Goat Anti-RERG Antibody. It is currently untested in blocking applications.

Comments

Peptide C-VKQAINKMLTKISS

Sequence

Shipping Ambient

Instructions

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade