

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05285 - IMMUNIZING PEPTIDE

SKU: EBP05285

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000953.2; NP 542158.1 ID

Peptide with sequence C-QDDGPAVERPSTEL, from the C Terminus of the protein sequence according to Immunogen

NP 000953.2; NP 542158.1.

The peptide was used in the production of EB05285 - Goat Anti-COX1 / PTGS1 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-QDDGPAVERPSTEL Sequence

Shipping

Grade

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png