

EBP05291-B - IMMUNIZING PEPTIDE

SKU: EBP05291-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_037390.2; NP_660183.1
Immunogen	Peptide with sequence CRESQSYLVEDLERS., from the C Terminus of the protein sequence according to NP_037390.2; NP_660183.1.
Product	The peptide was used in the production of EB05291-B - Goat Anti-ASC / TMS1, Biotinylated Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	CRESQSYLVEDLERS.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)