

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05294 - IMMUNIZING PEPTIDE

SKU: EBP05294

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_054728.2; NP_001107849.1 ID

Immunogen Peptide with sequence SQRPSRCSNPTPGP, from the C Terminus of the protein sequence according to NP_054728.2; NP_001107849.1.

The peptide was used in the production of EB05294 - Goat Anti-FOXP3 / SCURFIN Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

SQRPSRCSNPTPGP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade