

EBP05328 - IMMUNIZING PEPTIDE

SKU: EBP05328



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005477.2
Immunogen	Peptide with sequence C-IDAHQHKGAQNDGD, from the C Terminus of the protein sequence according to NP_005477.2.
Product Comments	The peptide was used in the production of EB05328 - Goat Anti-SRCASM / TOM1L1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-IDAHQHKGAQNDGD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)