

# EBP05340 - IMMUNIZING PEPTIDE

**SKU:** EBP05340



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_599011

**Immunogen** Peptide with sequence NSNEYCYVNSQQPS, from the C Terminus of the protein sequence according to NP\_599011.

**Product Comments** The peptide was used in the production of EB05340 - Goat Anti-Tim3 (mouse) Antibody. It is currently untested in blocking applications.

**Peptide Sequence** NSNEYCYVNSQQPS

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)