

# EBP05356 - IMMUNIZING PEPTIDE

**SKU:** EBP05356



---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                             |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.                                                                                                                                                                            |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                      |
| <b>Storage Instructions</b>  | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                                             |
| <b>Accession ID</b>          | NP_003921.2; NP_001290397.1                                                                                                                                                                                                                 |
| <b>Immunogen</b>             | Peptide with sequence C-GLVPKAYIMEMYDI, from the C Terminus of the protein sequence according to NP_003921.2; NP_001290397.1.                                                                                                               |
| <b>Product Comments</b>      | The peptide was used in the production of EB05356 - Goat Anti-SCAP2 / PRAP Antibody. It is currently untested in blocking applications.                                                                                                     |
| <b>Peptide Sequence</b>      | C-GLVPKAYIMEMYDI                                                                                                                                                                                                                            |
| <b>Shipping Instructions</b> | Ambient                                                                                                                                                                                                                                     |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)