

EBP05361 - IMMUNIZING PEPTIDE

SKU: EBP05361



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_037385.2; NP_056965.3; NP_001018067.1
Immunogen	Peptide with sequence C-AAAALLGGEEI, from the C Terminus of the protein sequence according to NP_037385.2; NP_056965.3; NP_001018067.1.
Product Comments	The peptide was used in the production of EB05361 - Goat Anti-DKK3 / REIC Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-AAAALLGGEEI
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)