

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05367 - IMMUNIZING PEPTIDE

SKU: EBP05367

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 054883.3; NP 659003.1; NP 001158077.1 ID

Peptide with sequence GKCRQLSRLCRNHCR, from the C Terminus of the protein sequence according to Immunogen

NP 054883.3; NP 659003.1; NP 001158077.1.

The peptide was used in the production of EB05367 - Goat Anti-Dysadherin Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide GKCRQLSRLCRNHCR Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade