

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05370 - IMMUNIZING PEPTIDE**

**SKU:** EBP05370

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** NP\_000894.1; NP\_001020604.1; NP\_001020605.1 ID

Peptide with sequence C-SIPTDNQIKARK, from the C Terminus of the protein sequence according to NP\_000894.1; Immunogen NP\_001020604.1; NP\_001020605.1.

The peptide was used in the production of EB05370 - Goat Anti-NQO1 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-SIPTDNQIKARK Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade