

EBP05393-B - IMMUNIZING PEPTIDE

SKU: EBP05393-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_570968.1; NP_001242909.1
Immunogen	Peptide with sequence CPSGRIRSVLHADEM., from the C Terminus of the protein sequence according to NP_570968.1; NP_001242909.1.
Product	The peptide was used in the production of EB05393-B - Goat Anti-AIBZIP / CREB3L4, Biotinylated Antibody. It is
Comments	currently untested in blocking applications.
Peptide Sequence	CPSGRIRSVLHADEM.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)