

EBP05396 - IMMUNIZING PEPTIDE

SKU: EBP05396



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1
Immunogen	Peptide with sequence C-APPPPTVSSLYL, from the C Terminus of the protein sequence according to NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1.
Product	The peptide was used in the production of EB05396 - Goat Anti-ADAM15 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-APPPPTVSSLYL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png