

EBP05399 - IMMUNIZING PEPTIDE

SKU: EBP05399



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_001100.3; NP_001157962.1 |
| Immunogen | Peptide with sequence C-QRKQGAGAPTAP, from the C Terminus of the protein sequence according to NP_001100.3; NP_001157962.1. |
| Product Comments | The peptide was used in the production of EB05399 - Goat Anti-MS2 / ADAM8 / CD156 Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-QRKQGAGAPTAP |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png |

DOCUMENTS

- [Data Sheet](#)