

EBP05401 - IMMUNIZING PEPTIDE

SKU: EBP05401



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006554.1
Immunogen	Peptide with sequence ATAETALPSISTLT-C, from the N Terminus of the protein sequence according to NP_006554.1.
Product Comments	The peptide was used in the production of EB05401 - Goat Anti-KLF1 / EKLF Antibody. It is currently untested in blocking applications.
Peptide Sequence	ATAETALPSISTLT-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)