

EBP05406 - IMMUNIZING PEPTIDE

SKU: EBP05406



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Shipped at ambient temperature, store at -20°C.

Instructions

Immunogen Peptide with sequence C-TAPPSHTGRHLW, from the C Terminus of the protein sequence according to .

Product The peptide was used in the production of EB05406 - Goat Anti-GPRC5B Antibody. It is currently untested in
Comments blocking applications.

Peptide
Sequence C-TAPPSHTGRHLW

Shipping
Instructions Ambient

Product
Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png